

Transitions™
XTRACTIVE®
FLASH TO MIRROR



BRING YOUR LENSES
TO THE NEXT LEVEL



Transitions™
XTRACTIVE®
FLASH TO MIRROR

BRING YOUR LENSES
TO THE NEXT LEVEL



Advantages
Vantaggi



Fashion colors

Colori fashion

A touch of fashion for all ages.
 Available in ten colors.
*Un tocco fashion per tutte le età.
 Disponibili in dieci colori.*

Protection


Protezione

The utmost comfort and protection.
Il massimo del comfort e della protezione.

Available colors

Colori disponibili

On Xtractive Brown base*
*Su base Xtractive Brown**



Bronze

Gold

Pink Gold

Red


On Xtractive Gray base*
*Su base Xtractive Gray**



Blue

Chrome

Green



Silver

Silver Green


Sky

* The images of the product are purely indicative. / * Le immagini del prodotto sono puramente indicative


FSV Organic Photochromic Index 1,5 Transitions® Xtractive® Flash to Mirror

Ranges / Gamme


Ø65 +/-



Ø70 +/-



Ø70 -/+



Available colors / Colori disponibili

Transitions® Xtractive® Flash to Mirror - **Blue**

Transitions® Xtractive® Flash to Mirror - **Bronze**

Transitions® Xtractive® Flash to Mirror - **Chrome**

Transitions® Xtractive® Flash to Mirror - **Gold**

Transitions® Xtractive® Flash to Mirror - **Green**

Transitions® Xtractive® Flash to Mirror - **Pink Gold**

Transitions® Xtractive® Flash to Mirror - **Red**

Transitions® Xtractive® Flash to Mirror - **Silver**

Transitions® Xtractive® Flash to Mirror - **Silver Green**

Transitions® Xtractive® Flash to Mirror - **Sky**

SFSV Organic Photochromic Index 1,5 Transitions® XTRActive® Flash to Mirror

Lenses Dimensions / Dimensioni Lenti

Base Curve (dpt)	True Curve (dpt)	Radius of Curvature	Back Curve	Diameter (mm)	C.T. (mm)	E.T. (mm)	SAG (30mm)	SAG (40mm)	SAG (50mm)	SAG (60mm)
0.50	0.44	1217.07.00	4.00	76	10.6	15.8	0.09	0.16	0.26	0.37
2.25	2.09	253.30.00	4.00	76	10.6	13.6	0.44	1.19	1.24	2.18
4.25	4.29	123.51.00	4.00	76	10.6	10.4	1.31	2.03	2.56	4.10
5.25	5.25	101.01.00	6.00	76	8.7	10.2	1.12	2.00	3.14	4.56
6.25	6.20	85.45.00	6.00	76	8.7	8.8	1.33	2.37	4.14	5.44
8.25	8.26	64.14.00	6.00	76	12.3	9.0	2.18	3.20	5.07	7.45

Available colors / Colori disponibili

Transitions® Xtractive® Flash to Mirror - **Blue**

Transitions® Xtractive® Flash to Mirror - **Bronze**

Transitions® Xtractive® Flash to Mirror - **Chrome**

Transitions® Xtractive® Flash to Mirror - **Gold**

Transitions® Xtractive® Flash to Mirror - **Green**

Transitions® Xtractive® Flash to Mirror - **Pink Gold**

Transitions® Xtractive® Flash to Mirror - **Red**

Transitions® Xtractive® Flash to Mirror - **Silver**

Transitions® Xtractive® Flash to Mirror - **Silver Green**

Transitions® Xtractive® Flash to Mirror - **Sky**

Organic - Single Vision

Index 1,5 Transitions® XTRActive®

Flash to Mirror



Ranges / Gamme

-10.00	-10.00	-10.00	-10.00	-5.00
∅ 50	∅ 55	∅ 60	∅ 65	∅ 70
+6.00	+6.00	+6.00	+6.00	+6.00
CYL UP TO +6.00				CYL UP TO +4.00


Available colors / Colori disponibili

- Transitions® Xtractive® Flash to Mirror - **Blue**
- Transitions® Xtractive® Flash to Mirror - **Bronze**
- Transitions® Xtractive® Flash to Mirror - **Chrome**
- Transitions® Xtractive® Flash to Mirror - **Gold**
- Transitions® Xtractive® Flash to Mirror - **Green**
- Transitions® Xtractive® Flash to Mirror - **Pink Gold**
- Transitions® Xtractive® Flash to Mirror - **Red**
- Transitions® Xtractive® Flash to Mirror - **Silver**
- Transitions® Xtractive® Flash to Mirror - **Silver Green**
- Transitions® Xtractive® Flash to Mirror - **Sky**

Organic - Progressive Index 1,5 Transitions® XTRActive® Flash to Mirror



Ranges / Gamme



Available colors / Colori disponibili

- Transitions® Xtractive® Flash to Mirror - **Blue**
- Transitions® Xtractive® Flash to Mirror - **Bronze**
- Transitions® Xtractive® Flash to Mirror - **Chrome**
- Transitions® Xtractive® Flash to Mirror - **Gold**
- Transitions® Xtractive® Flash to Mirror - **Green**
- Transitions® Xtractive® Flash to Mirror - **Pink Gold**
- Transitions® Xtractive® Flash to Mirror - **Red**
- Transitions® Xtractive® Flash to Mirror - **Silver**
- Transitions® Xtractive® Flash to Mirror - **Silver Green**
- Transitions® Xtractive® Flash to Mirror - **Sky**



Customer assistance

Customer service e-mail
customer.assistance@lenses.it

Opening times

From monday to friday from 09:00 to 18:00

Assistenza clienti

e-mail assistenza clienti
assistenza.clienti@lenses.it

Orari

Dal lunedì al venerdì dalle 09:00 alle 18:00

www.lenses.it